

# KYLT<sup>®</sup>

**SAN**<sup>®</sup>  
VET



## **KYLT<sup>®</sup> SALMONELLA HADAR**

**REAL-TIME PCR DETECTION**

**i For in vitro Veterinary Diagnostics only.**

[WWW.KYLT.EU](http://WWW.KYLT.EU) | [WWW.SAN-VET.COM](http://WWW.SAN-VET.COM)

**Passion for Innovation**



## DIRECTION FOR USE

qPCR.Salm.Hadar.02, Rev002, November 2024

# KYLT<sup>®</sup>

## SALMONELLA HADAR

### Real-Time PCR Detection

#### 1. GENERAL





Characteristic	Description
Article name	Kylt <sup>®</sup> S. Hadar LD 100, Kylt <sup>®</sup> S. Hadar LD 25
Name of pathogen	Salmonella Hadar
Caused disease	Salmonellosis
Organism	Bacteria
Target molecule	DNA
Technology	Real-Time PCR
Host-group	Poultry (Chicken, Turkey)
Sample prerequisite	Salmonella spp. 2.0
Sample types	<ul style="list-style-type: none"><li>pre-enrichment samples positively tested for Salmonella spp. (Ct &lt; 35), Salmonella spp. positive colony material from cultural processes (e.g. DIN EN ISO 6579)</li></ul>
Target FAM channel (520 nm)	Salmonella group O:8 (C2-C3)
Target Cy5 channel (670 nm)	Salmonella H-phase 1 (z10)
Target TXR channel (620 nm)	Salmonella H-phase 2 (e, n ,x)
Temperature profile	Kylt <sup>®</sup> Profile I or Kylt <sup>®</sup> Profile II

These kits were developed for use by trained laboratory personnel following standardized procedures. This Direction For Use must be followed strictly.

## 2. REAL-TIME PCR

During the Real-Time PCR the target genes are amplified by respective primer pairs in the Polymerase Chain Reaction (PCR). Amplified target sequences are detected via fluorescently labeled probes during the PCR reaction in real-time (Real-Time PCR). Their emitted fluorescence is separately optically measured by the Real-Time PCR thermal cycler. The pathogen-specific status of a sample can be evaluated by considering all amplified targets per sample and the negative and positive controls per run.

## 3. KIT CONTENTS

Reagent	Lid color	100 Reactions Article no. 31547	25 Reactions Article no. 31548	Storage temperature
2x qPCR-MIX	 Transparent	4× 280 µl	1× 280 µl	<= -18°C
Primer-Probe-Mix	 Green	4× 450 µl	1× 450 µl	<= -18°C
Positive Control	 Red	4× Lyoph. a final 50 µl	2× Lyoph. a final 50 µl	<= -18°C
Negative Control	 Blue	1× 1 ml	1× 1 ml	<= -18°C

## 4. STORAGE REQUIREMENTS

- If necessary, prepare aliquots of the reagents to limit freeze-thaw cycles to 3.
- The components are to be used within the indicated shelf life (see box label).
- Components from different batches may not be mixed or interchanged.
- The **Primer-Probe-Mix** needs to be stored and handled protected from light.

## 5. REAGENT PREPARATION

- Before its first use, rehydrate each vial of **Positive Control**: add 50 µl of Negative Control per vial, briefly incubate at room temperature and mix thoroughly by repeated vortexing.
- Before first use, rehydrate the lyophilized **Primer-Probe-Mix**: add 150 µl of the Negative Control per vial, briefly incubate at room temperature and mix by pulse-vortexing.

## 6. NECESSARY EQUIPMENT, DEVICES AND CONSUMABLES

- Real-time PCR cycler capable of detecting the appropriate fluorescence wavelength. (Note that the default normalization option against ROX (e.g. using ABI cyclers) must be disabled).
- Compatible PCR-plates strips or single tubes.
- Purification Kit yielding sufficiently high concentration of inhibitor-free DNA/RNA (e.g. Kylt RNA / DNA Purification products)."
- Table top microcentrifuge.
- Vortex mixer.
- Adjustable Micropipettes covering the appropriate volume range.
- Matching PCR-clean pipette tips with filters.
- Certified Nuclease-free (PCR-clean) consumables.
- Powder-free gloves to be worn during the entire setup and changed in case of contamination.

## 7. CONTROL REACTIONS

- Each PCR run must include a **Positive Control** to monitor the specificity and efficiency of the reagents and the reaction itself, including the performance of the Real-Time PCR and of the Real-Time PCR thermal cycler.
- Each PCR run must include a **Negative Control** to ensure the absence of contaminations.

## 8. REACTION SETUP

- Before each use, briefly vortex and spin down the used reagents.
- The total number of reactions is the number of samples plus one Positive Control and one Negative Control per run. In case of combination this setup with a Kylt Standard, please add appropriate numbers of reactions.
- Prepare the Master-Mix using the components listed below.
- Vortex, spin down and add 16 µl of the finalized Master-Mix to each of the PCR tubes or plate wells ("cavities").

Reagent	Volume (µl)	
	per Reaction	e.g. n= 7
2x qPCR-Mix	10 µl	70 µl
Primer-Probe-Mix	6 µl	42 µl
<b>Total Volume Master-Mix</b>	<b>16 µl</b>	<b>112 µl</b> (dispense 16 µl per reaction)
DNA (Negative Control / sample DNA / Positive Control)	4µl	
<b>Total Volume per Reaction</b>	<b>20 µl</b>	

- Keep exposure of the 2x RT-qPCR-Mix, Primer-Probe-Mix and prepared Master-Mix to (sun)light as short as possible and return it back to appropriate storage temperature right after application.
- Avoid the formation of bubbles when pipetting the mix, samples and controls.
- Add 4 µl of the **Negative Control** to the corresponding cavity and seal it individually, if possible.
- Add 4 µl of each **sample** to the corresponding cavities and seal them individually, if possible.
- Add 4 µl of the **Positive Control** to the corresponding cavity.
- Once all reactions are set up, seal the cavities and briefly spin down.
- Place the cavities in the Real-Time PCR thermal cycler and run the test with either Kylt Profile II or Kylt Profile I.
- With Kylt Profile II this and most other Kylt qPCR detection methods can be carried out simultaneously in a single PCR run.
- Kylt Profile I allows for combined run of this and most other Kylt RT-qPCR detection methods as well as Kylt PCR detection products.
- In the event of a combined Real-Time (RT-)PCR run, make sure all necessary channels are detected.
- Please follow the specified instructions of your Real-Time PCR thermal cycler as recommended by the manufacturer.

### Kylt® Profil II

Step No	Description	Temperature	Duration
1	Activation of Polymerase	95 °C	10 min
2	Denaturation	95 °C	15 sec
3	Annealing & Extension	60 °C	1 min
4	Fluorescence Detection	channels / wave lengths see page 2	

} 42 cycles

### Kylt® Profil I

Step No	Description	Temperature	Duration
1	Reverse Transcription	50 °C	10 min
2	Activation of Polymerase	95 °C	1 min
3	Denaturation	95 °C	10 sec
4	Annealing & Extension	60 °C	1 min
5	Fluorescence Detection	channels / wave lengths see page 2	

} 42 cycles

## 9. GENERAL PCR-RUN EVALUATION AND QUALITY CONTROL

- An ideal positive PCR curve starts with a linear phase which then increases exponentially and tapers off to a plateau phase.
- The baseline is the mean value of a large section of the beginning of the experiment which is subtracted from all fluorescence values to show only the net increase of fluorescence per well and channel.
- Automated evaluation of the respective Cycler Software may be used. Please take care to identify possible artifacts.
- The threshold, if adjusted manually, should be set close enough to the baseline to include all curves that show a clear exponential phase but to exclude all unspecific fluorescence increases.
- The intersection between the curve and the threshold is the Ct-value. The lower the Ct-value, the higher the concentration of the target molecule in the sample at the beginning of the test.

## 10. TEST EVALUATION - CONTROL REACTIONS

- Only Salmonella spp. positive samples (Ct-value < 35) or isolates can be analyzed with this PCR kit.
- The Real-Time PCR run is only valid if the curves of the control reactions can be evaluated as follows:

Control Reactions	Pathogen-Specific Channels					
	FAM		Cy5		TXR	
Negative Control	negative	Ct > 42	negative	Ct > 42	negative	Ct > 35
Positive Control	<b>positive</b>	<b>Ct &gt; 15 und ≤ 35</b>	<b>positive</b>	<b>Ct &gt; 15 und ≤ 35</b>	<b>positive</b>	<b>Ct &gt; 15 und ≤ 35</b>

## 11. TEST EVALUATION - SAMPLES

Target	Channel	Signal				
<i>Salmonella</i> group O:8 (C2-C3)	FAM	<b>positive</b>	<b>positive</b>	negative	negative	negative
<i>Salmonella</i> H-phase 1 (z10)	Cy5	<b>positive</b>	negative	<b>positive</b>	negative	negative
<i>Salmonella</i> H-phase 2 (e, n ,x)	TXR	<b>positive</b>	negative	negative	<b>positive</b>	negative
The sample is <i>Salmonella Hadar</i>		<b>positive</b>	negative	negative	negative	negative

- The **pathogen-specific channel** is positive: Ct-value  $\leq$  42. **Exception:** The **TXR channel** is positive: Ct-value  $<$  35.
- As double positive samples (samples containing more than one *Salmonella* serovar) can not be excluded it is strongly recommended to confirm PCR positive results by cultural, immunological or biochemical methods according to your local regulations.
- **Recommendation:** In the case of an inhibited sample the test may be repeated with a dilution of the DNA preparation at e.g. 1:10. The Negative Control is used as the diluting agent. Preferably, the entire DNA preparation process is repeated using Kylt RNA/DNA Purification products or appropriate alternative.
- Convenient and reliable sample data entry, Real-Time PCR start, final qualitative analysis and documentation can be conducted with the Kylt Software, please inquire.

## 12. ORDERING INFORMATION

For a fast and efficient service please send your order to [orders.kylt-de@san-group.com](mailto:orders.kylt-de@san-group.com) and please provide the following information:

- Delivery address
- Invoice address
- Purchaser contact telephone number
- End user name and telephone number (if different)
- Purchase order number
- Product name and catalogue number
- Quantity and size of products
- Indicate if your account is VAT exempt

## 13. REVISION HISTORY

Revision	Status	Amendments
Rev002	November 2024	changed kit content (2x qPCR-Mix and Primer-Probe-Mix instead of Reaction-Mix)

**Production:** SAN Group Biotech Germany GmbH · Muehlenstrasse 13 · 49685 Hoeltinghausen · Germany  
[www.kylt.eu](http://www.kylt.eu) · [kylt-de@san-group.com](mailto:kylt-de@san-group.com)

The development, manufacturing and distribution of Kylt® In-Vitro Diagnostics is certified to ISO 9001:2015.

Kylt® is a registered trademark.

For veterinary use only. For in vitro use only. Regulatory requirements vary by country, not all of the products described herein may be available in your geographic area.

© 2024 SAN Group Biotech Germany GmbH. The trademark mentioned herein is the property of SAN Group Biotech Germany GmbH or their respective owners.






# PROTOCOL AT A GLANCE

## REAL-TIME PCR SETUP


1

Prepare a Master-Mix




2

Pulse-vortex and spin down



3

Dispense Master-Mix and add 4 µl NC, sample DNA/RNA, PC



4

Seal cavities, spin down (recommended), and start cycler



5

Analysis

